



Foetal Bovine Serum

- ◆ **Australian Origin FBS**
- ◆ **European Directorate of Quality Medicines (EDQM) Certified**
- ◆ **Produced by Bovogen Biologicals in their fully validated 'state of the art' filtration facility**
- ◆ **Superior quality at competitive price**
- ◆ **FBS source fully traceable and documented**

Bovogen Biologicals is a leading manufacturer and supplier of genuine Australian foetal bovine serum (FBS) to the global life sciences community.

Australian origin Foetal Bovine Serum is highly sought after due to Australia's freedom from all major epidemic diseases of livestock, including Foot and Mouth Disease.

This status is collectively administered by The World Organisation for Animal Health (OIE), the Food and Agriculture Organisation of the United Nations (FAO) and the World Health Organisation (WHO), who are responsible for monitoring the international epidemiology of livestock diseases.

Bovogen Biologicals' Foetal Bovine Serum is manufactured in our purpose built sterile filtration facility, under strict GMP conditions, which ensures that it meets the highest quality standards.



Foetal Bovine Serum

All foetal bovine serum is collected from registered export abattoirs under the supervision of Government Veterinary authorities. In line with our quality systems, the origin of Bovogen FBS is fully traceable back to the abattoir and is supported by the availability of full collection and manufacturing records, leaving no doubt as to its' source.

The production process is closely monitored from blood collection to sterile filtration in our fully validated and dedicated state-of-the-art filtration facility. Notably this process includes a closed loop system to ensure that the final product has low endotoxin and haemoglobin.

In further support of our quality manufacturing practices, Bovogen Biologicals received a Certificate of Suitability from the European Directorate of Quality Medicines for its range of Australian origin foetal bovine serum (FBS) products in June, 2008.

This certification meets the criteria described in the current version of the monograph "Products with Risk of Transmitting Agents of Animal Spongiform Encephalopathies No. 1483 of the European Pharmacopoeia.

This Certificate of Suitability now enables Bovogen Biologicals' foetal bovine serum to be used for the manufacture of medicines in the European Union.

As an Australian manufacturer and supplier of Australian origin Foetal Bovine Serum, Bovogen Biologicals' mission is to guarantee that the quality of our products will meet and exceed our specifications each and every time.

We take pride in our ability to surpass customer expectations for product quality, delivery and service and we are committed to ensuring that Bovogen Biologicals is not just a supplier but an industry partner for our customers.

SPECIFICATIONS

Product	Foetal Bovine Serum
Catalogue No	SFBS
Size	500ml
Animal Origin	Australian
Format	Sterile Filtered by 3 x 0.1µm filters which conform to FDA regulations 211.72.
Reported Parameters	pH, Osmolality, Haemoglobin, Total Protein, Total IgG, Endotoxin, Sterility, Mycoplasma, Biochemistry
Virus Isolation Test	BVD, IBR, PI3
Certification	European Directorate of Quality Medicine's (EDQM) Certificate of Suitability

To request your free Foetal Bovine Serum sample, please complete the information below and return to Bovogen Biologicals on fax no: +613 9336 3644 or email: info@bovogen.com

Name: _____

Title: _____

Organization: _____

Organization Address: _____

Building/Room/Laboratory: _____

City: _____ State: _____

Country _____ Zip/Postcode: _____

Business Telephone: _____ Mobile: _____

Email Address: _____



Bovogen Biologicals Pty Ltd

Head Office

71 Ogilvie Street
Essendon, Victoria 3040
Australia
Phone +613 9337 6954
Fax +613 9337 0738
info@bovogen.com

Manufacturing Facility

12 Williams Avenue
East Keilor, Victoria 3033
Australia
Phone +613 9336 3622
Fax +613 9336 3644
www.bovogen.com

JULY 2010



Bovogen Biologicals is a certified member of the Australia made campaign.